

RECEIVED  
CENTRAL FAX CENTER

DEC 12 2005

**IN THE UNITED STATES PATENT AND TRADEMARK OFFICE**

Application of E. C. Leuthardt  
Serial No.: 10/735,474  
Confirmation No. 2473  
Filed: December 12, 2003  
Entitled: BRAIN COMPUTER INTERFACE

Group Art Unit 3736  
Examiner: P. Mallari

**DECLARATION UNDER 37 C.F.R. § 1.131**

I, Eric C. Leuthardt, declare and state as follows:

1. I am a co-inventor of the subject matter claimed in U.S. Patent Application 10/735,474 entitled, "BRAIN COMPUTER INTERFACE," which was filed December 12, 2003.
2. U.S. Patent Application Publication No. US 2005/0113744 A1 was cited against claims of the subject application in an Office Action dated July 11, 2005. U.S. Patent Application Publication No. US 2005/0113744 A1 has a priority date of November 21, 2003, which is twenty days before the priority date of the subject application.
3. Attached as Exhibit A is a copy of my laboratory notes dated March 7, 2003 showing results for a Screening Protocol conducted with subject JJ using the Brain Computer Interface of the subject patent application.
4. Attached as Exhibit B is a copy of my laboratory notes dated May 21, 2003 showing results for a Screening Protocol conducted with subject DW using the Brain Computer Interface of the subject patent application.
5. Based on the attached exhibits, it is evident I and my co-inventors constructively reduced to practice the invention recited in the claims of the subject application before the priority date of U.S. Patent Application Publication No. US 2005/0113744 A1 and diligently pursued filing U.S. Patent Application 10/735,474 during a period of time beginning before the priority date of U.S. Patent Application Publication No. US 2005/0113744 A1 and lasting at least until the filing date of Application Ser. No. U.S. Patent Application 10/735,474.

4. I hereby declare that the statements made of my own knowledge are true and that all statements made on information or belief are believed to be true. I acknowledge that willful false statements and the like are punishable by fine or imprisonment or both (18 U.S.C. § 1001) and may jeopardize the validity of this application or any patent issuing therefrom.

Date:

11/10/05  
Eric C. Leuthardt, M.D.

## EXHIBIT A

RECEIVED  
CENTRAL FAX CENTER

DEC 12 2005

Purpose:

Screening Protocol St. Louis 03/07/2003

Initials: JJ

Date: 3/7/02Technician: Eric Leubandt

Run	Duration	Task	.prm File	Comments
①	1 min	resting with eyes open	base60	
②	1 min	resting with eyes closed	base60	noised @ arm
③	2 min	opening and closing left/right hand	horizontal	Did well
④	2 min	imagining opening and closing left/right hand	horizontal	Focused well
⑤	2 min	protrude tongue in and out	onetarget	Protruded tongue
⑥	2 min	imagining to protrude tongue in and out	onetarget	
⑦	2 min	saying a word "move"	onetarget	Seemed more distracted
⑧	2 min	imagine saying a word "move"	onetarget	Said word only once
⑨	3 min	joystick task	joystick	
⑩	2 min	opening and closing left/right hand	horizontal	5: no excessive arm movement
⑪	2 min	imagining opening and closing left/right hand	horizontal	
⑫	2 min	protrude tongue in and out	onetarget	inconsistent - did it consistently
⑬	2 min	imagining to protrude tongue in and out	onetarget	Seemed hyperkinetic
⑭	2 min	saying a word "move"	onetarget	Seemed inconsistent
⑮	2 min	imagine saying a word "move"	onetarget	Seemed focused
⑯	3 min	joystick task	joystick	did well
⑰	2 min	opening and closing left/right hand	horizontal	inconsistent

18	2 min	Imagining opening and closing left/right hand	horizontal	
19	2 min	protrude tongue in and out	onetarget	
20	2 min	imagining to protrude tongue in and out	onetarget	protrudes
21	2 min	saying a word "move"	onetarget	
22	2 min	Imagine saying a word "move"	onetarget	
23	3 min	joystick task	joystick	

VERY PERSEVERATIVE → goes say 1 step points

- UL:
- JJ is finally recovered enough to complete all the tasks required for the screening protocol. The protocol is attached above.
  - He is somewhat tired and pained here with the tasks.
  - Data acquired is good.
  - We will attempt to do long tomorrow.
  - Note: we have cancelled building the computer useful pin out diagrams for simulation in back

RE: THE FILTERS IN THE ACQUISITION COMPUTER AND THE ANALOG COMPUTER WERE ON. THIS PROBABLY MEANS THAT HIGHER FREQUENCIES ABOVE 70 Hz WERE FILTERED OUT.

## EXHIBIT B

Session DW 001

Date: 5/21/03 - 5/22

Screening Session: 001

Run	Task	.prm file	Comments
001 ✓	Resting eyes open	Base60	did well
002 ✓	Resting eyes closed	Base60	did well
003 ✓	Open closing right/left hand	horizontal	did well
004 ✓	Imagine Open closing right/left hand	horizontal	too fast
005 ✓	Protrude tongue in and out	Onetarget	did well
006 ✓	Imagine protrude tongue in and out	Onetarget	did well
007 ✓	Say word move	Onetarget	did well
008 (5/22/03)	Imagine saying word move "us"	Onetarget	Network crashed (5/22)
009 ✓	joystick	joystick	did well
010 ✓	Open closing right/left hand	horizontal	did well
011 ✓	Imagine Open closing right/left hand	horizontal	did well
012 ✓	Protrude tongue in and out	Onetarget	did well
013 ✓	Imagine protrude tongue in and out	Onetarget	did well
014 ✓	Say word move	Onetarget	did well
015 ✓	Imagine saying word move	Onetarget	did well
016 ✓	joystick	joystick	did well
017 ✓	Open closing right/left hand	horizontal	did well
018 ✓	Imagine Open closing right/left hand	horizontal	did well
019 ✓	Protrude tongue in and out	Onetarget	did well
020 ✓	Imagine protrude tongue in and out	Onetarget	telling asump
021 ✓	Say word move	Onetarget	did well
022 ✓	Imagine saying word move	Onetarget	did well
023 ✓	joystick	joystick	did well

024

joystick

did well

did well

## Pre-Testing Check List — 5/21/05 : DW

X	Unplug Bed
X	Unplug Fan and turn off TV
X	Use cortical lead as reference (most inactive lead) requires tech assistance
X	Change the HIGH filter settings on the acquisition computer (in EEG monitor lab) to OFF
-	Make sure that the acquisition frequency is set at 500 HZ
	Uncheck the analog printer filter settings on the BCI computer

Names of Individuals involved :

Eric Leuthardt

Gerwin Schalk

Dan Moran

Note = reference = scalp

## Post Testing Check List

- |  |   |
|--|---|
|  | Make sure name and run are correct following each run                 |
|  | Make sure patient is reconnected to local network and being monitored |

Note : In checking parameter files we noted that the last one was sampled @ 500hz. Therefore we need to rev data accordingly to 300hz. Because reanalysis analyzed the data according to

etc. add to check list - the reversed pm parameters

## SOURCE

- \* Sampling rate 500
- Software channels 32
- Sample Block size 20

## MEM FACTOR

Bandwidth 2

Model order 18

MEM window 7

Channel 1-18

**This Page is Inserted by IFW Indexing and Scanning  
Operations and is not part of the Official Record**

**BEST AVAILABLE IMAGES**

Defective images within this document are accurate representations of the original documents submitted by the applicant.

Defects in the images include but are not limited to the items checked:

- ☐ **BLACK BORDERS**
- ☐ **IMAGE CUT OFF AT TOP, BOTTOM OR SIDES**
- ☐ **FADED TEXT OR DRAWING**
- ☐ **BLURRED OR ILLEGIBLE TEXT OR DRAWING**
- ☐ **SKEWED/SLANTED IMAGES**
- ☐ **COLOR OR BLACK AND WHITE PHOTOGRAPHS**
- ☐ **GRAY SCALE DOCUMENTS**
- ☐ **LINES OR MARKS ON ORIGINAL DOCUMENT**
- ☐ **REFERENCE(S) OR EXHIBIT(S) SUBMITTED ARE POOR QUALITY**
- ☐ **OTHER:** \_\_\_\_\_

**IMAGES ARE BEST AVAILABLE COPY.**

**As rescanning these documents will not correct the image problems checked, please do not report these problems to the IFW Image Problem Mailbox.**